

Microstate EEGLAB toolbox

Andreas Pedroni, Andreas Poulsen, Nicolas Langer, Lars Kai
Hansen

Methods of Plasticity Research, University of Zurich
Department of Applied Mathematics and Computer
Science, Technical University of Denmark



**University of
Zurich**^{UZH}





Technical University of Denmark

Classical EEG analysis

EEG

- Time (ERP)
- Frequency (FFT)
- Time-Frequency (Wavelet)



Typically many electrodes

Which electrode(s) should I analyze?

→ Problem I: multiple testing, cherry picking,

→ Problem II: missing out on information



A state-dynamics view of EEG signals

One approach for analyzing EEG-signals in a more holistic fashion is to view it as **states** with **dynamics**


- **States** = Description of the EEG-signal at any time t (complexity, synchronicity, etc.)
- **Dynamics** = How do the states behave over time (duration, frequency of occurrence, sequence)

A topographical view on EEG

Coherent firing of a large number of neurons produces a defined voltage field on the scalp (=topography)




Electrodes



Time

A topographical view on EEG




- **Change of the topography** indicates change in activation of networks
- **Change of the amplitude** indicates change in the number of coherently active neurons

Topographical analysis describes the **dynamics** of the functional **states** of neuronal networks

The birth of EEG microstates: EEG topographies are „quasi-stable“ over time






Data: 1500ms, topographies at GFP peaks



Observation:


80 – 150ms, quasi-stable topographical configurations



„The spatially stationary segments might be basic building blocks of brain information processing, possibly operationalizing consciousness time and offering a common phenomenology for spontaneous activity and event-related potentials“ (Lehmann et al., 1987)

EEG Microstates in spontaneous data

- During resting (closed eyes):
- Typically 4 microstate class topographies (explain ~80 % of the topographical variance in the EEG)
- Microstates last on average between 80 – 120ms




The basic steps of a EEG microstates Analysis

1. Preprocessing
2. **Microstates segmentation** (*pop_micro_segment*)
Goal: Find the prototypical topographies
(=microstate-classes, prototype maps)
3. **Back-fitting** (*pop_micro_fit*)
Goal: Label the data of interest with respect to the best fitting Microstate-class
4. **Microstates statistics** (*pop_micro_stats*)
Goal: Calculate for the EEG segments of a respective MS-class the parameters of interest
5. Inferential statistical analysis of parameters

Microstates analysis in spontaneous vs. event-related data

- **Spontaneous data (e.g. resting data):**
- Ignore polarity of topography
- Modified K-means (Pascual-Marqui et al., 1995)
- Global Field Power (GFP) peaks
 - topographies with high signal-to-noise ratio
 - topographies tend to be stable around GFP peaks
- **Event related data:**
- Account for polarity of topography
- K-means
- Use (grand)averaged data for segmentation



Microstates segmentation and back-fitting


Segmentation:

The segmentation aims to find the most prototypical topographies (=microstate-classes, prototype maps) of the data

- use GFP peaks
- Modified k-means clustering
- How many clusters? Cross validation, other criteria...

Back-fitting:

Assign the microstate-class with highest spatial correlation to each timeframe of EEG (with minimum segment duration)



without min duration:



The Output: EEG microstates statistics

- Average duration
- Mean number of occurrences per second
- Percentage time covered by each state
- Global explained variance of each state
- Average GFP during a state
- Transition probabilities between states

Functional significance of EEG microstates

Emergence of microstates:


Improbable: A microstate is produced by a small, local group of neurons that are transiently active together (Only small number of MS, well defined temporal structure)

Probable: A microstate is produced by a large network of activity, spanning over the large parts of the cortex

Connection to fMRI resting state networks:

- When the microstate time series is convoluted with the resting-state fMRI BOLD signal, individual **microstate maps correlate with the activity of particular resting state networks** (Britz et al., 2010; Musso et al., 2010; Yuan et al., 2012).
- MS time series are linked to fMRI signal via scale free dynamics, suggesting that resting-state brain activity at the **EEG and fMRI timescales reflect the same underlying neurophysiological processes** (Van de Ville et al., 2010)

Now to the practice...

- Resting EEG data of children
- Filtered from 2 to 20 Hz
- Preprocessed with Automagic 
- Plugin: <https://github.com/atpoulsen/Microstate-EEGLab-toolbox>
- <https://github.com/DynAge/brainhack-zh/blob/master/misc/MicrostateBrainHackSingle.m>
- Data: the same as in Automagic.



Further reading...

Topographic ERP Analyses: A Step-by-Step Tutorial Review by Murray et al., 2008

Microstates in resting-state EEG: Current status and future directions by Khanna et al., 2015

Electrical Neuroimaging by Michel et al., 2009

Email: andreas.pedroni@uzh.ch